




Oncomine Focus Assay

Single Panel for Testing 18 Types of Solid Tumors

Oncomine focus

- Hot spot screening of relevant, druggable mutations, copy number changes and fusions for various Cancers
- Analysing of 1000 variants across 52 genes for --cancers
- ACMG/AMP guidelines-based reporting with drugs and target information
- NGS based assay eliminate need of multiple technology FISH, IHC, Realtime
- High Sensitivity and Specificity
- Identifies druggable targets in various cancers and unknown primary

Workflow



Actionable Report

Oncomine Knowledge base
Reporter

List of genes

SNV's; Indels	<p><i>AKT1, ALK, AR, BRAF, CDK4, CTNNB1, DDR2, EGFR, ERBB2, ERBB3, ERBB4, ESR1, FGFR2, FGFR3, GNA11, GNAQ, HRAS, IDH1, IDH2, JAK1, JAK2, JAK3, KIT, KRAS, MAP2K1, MAP2K2, MET, MTOR, NRAS, PDGFRA, PIK3CA, RAF1, RET, ROS1, SMO</i></p>
CNV's	<p><i>AKT1, ALK, AR, BRAF, CCND1, CDK4, CDK6, EGFR, ERBB2, FGFR1, FGFR2, FGFR3, FGFR4, KIT, KRAS, MET, MYC, MYCN, PDGFRA, PIK3CA</i></p>
Fusions	<p><i>ABL1, AKT3, ALK, AXL, BRAF, EGFR, ERBB2, ERG, ETV1, ETV4, ETV5, FGFR1, FGFR2, FGFR3, MET, NTRK1, NTRK2, NTRK3, PDGFRA, PPARG, RAF1, RET, ROS1</i></p>



Thyroid Cancer 18

AKT, BRAF, CTNB1, HRAS, IDH1, KRAS, NRAS, NTRK1, RET

Testicular Carcinoma 17

BRAF, CTNB1, FGFR3, HRAS, KIT, KRAS, NRAS

Soft Tissue Sarcoma 16

CCND1, CDK4, ERG, FGR1, HRAS, KRAS, NRAS, PIK3CA

Skin Basal Cell Carcinoma 15

CTNB1, HRAS, KRAS, NRAS, PIK3CA

Prostate Cancer 14

AR, BRAF, ERG, ETV1, ETV4, ETV5, FGFR1, HRAS, IDH1, KRAS, MYC, RAF1

Pancreatic Cancer 13

CDK6, CTNB1, FGFR1, KRAS, MYC

Ovarian Cancer 12

BRAF, CCND1, CTNB1, FGFR1, KIT, KRAS, MYCN, PIK3CA


Melanoma Cancer 11

BRAF, CCND1, CDK4, GNA11, GNAQ, IDH1, KRAS, MAPK21, MYC, NRAS

Liver Cancer 10

BRAF, CCND1, CTNB1, IDH1, IDH2, KRAS, MYC, NRAS

List
key dru
and d
mutation
differen



01 Bladder Cancer

AKT1, CCND1, EGFR,
FGFR1, FGFR3, HRAS,
KRAS, MUCM PIK3CA



02 Breast Cancer

AKT1, CCND1, ERBB2,
FGFR1, MYC, NTRK3,
PIK3CA



03 Colon Cancer

KRAS, NRAS, HRAS,
BRAF, PIK3CA, AKT1



04 Endometrial Cancer

AKT1, CCND1, CTNB1,
ERBB2, ESR1, FGFR1, FGFR2,
KRAS, MYC, NRAS, PIK3CA



05 Esophageal Cancer

CCND1, CDK6, CTNB1,
EGFR2, FGFR1, FGFR2, KRAS,
MAP2K1, MYC, PIK3CA



06 Gastric Cancer

AR, BRAF, CCND1, CDK4, CDK6,
EGFR, ERBB2, ERBB3, FGFR3,
KRAS, MET MYC, PIK3CA



07 GIST Cancer

BRAF, CCND1, KIT,
PDGFRA, PIK3CA



09 Lung Cancer

EGFR, BRAF, MET,
RET, ERBB2, KRAS,
ALK, RET, ROS, NTRK




08 Head & Neck Cancer

CCND1, CDK6, EGFR,
ERBB2, FGFR1, FGFR3,
HRAS, KRAS, MET,
NRAS, PIK3CA




SNV and CNV profile of oncomine focus genes in Pan-cancer data



All types of solid tumors showing the presence of mutation and copy number variation, in the genes present in oncomine focus assay. Hence, this assay best suited to identify druggable mutation, driver mutation or CNV changes and prognostic bio-markers.

Analyzed using: cbiportal.org.

Oncomine Fusion assay genes profile in Pan-cancer data



Multiple tumors such as Bone sarcoma, breast cancer, CRC, cholangiocarcinoma, gastric cancer, Head and neck cancer, melanoma, lung cancer, prostate cancer, salivary gland carcinoma, thyroid cancer, uterine sarcoma are showing fusion drivers in >2% of the cases, that can be detected by Oncomine focus assay.

Analyzed using: cbiportal.org.

Evidence based variant classification by Oncomine Reporter

Tier I (Therapeutic, Diagnostic, Prognostic)

Variants of strong clinical significance

Level A

FDA Approved

Professional Guidelines

Level B

Well powered studies

Tier II (Therapeutic, Diagnostic, Prognostic)

Variants of potential clinical significance

Level C

FDA Approved
In different tumor,
small published
studies

Level D

Preclinical trials


Tier III (Variants of unknown clinical significance)

Not observed at a significant allele frequency in the general or specific subpopulation databases, or pan-cancer or tumor-specific variant databases

Tier IV

Benign or likely benign variants

Assay sensitivity and specificity



- Can identify mutation with 5% variant allele frequency (5% mutant DNA, among 95% normal DNA).
- LOD for RNA fusion is 255 copies of RNA.
- RNA fusions and CNV changes detected by NGS has equal sensitivity to FISH. Hence multiple testing is not required.

List of druggable targets

List of FDA approved drugs covered

Biomarker	Drugs	Disease
BRAF V600E; BRAF V600K	Vemurafenib, Cobimetinib, Dabrafenib, Trametinib, Atezolizumab, Encorafenib	Melanoma, NSCLC, CRC
ALK Fusion, EML4-ALK	Crizotinib, Ceritinib, Alectinib, Brigatinib, Lorlatinib	NSCLC
c-KIT	Imatinib, Sunitinib, Regorafenib, Ripertinib	GIST
EGFR	Cetuximab, Mobocertinib, Afatinib, Gefitinib, Amivantamb, Osimertinib, Erlotinib	Lung cancer
EGFR	Panitumumab, Dacomitinib	CRC
ERBB2	Fam-trastuzumab Deruxtecan-nxki	NSCLC
ERBB2/HER2 amplification	Trastuzumab, Pertuzumab, Ado-trastuzumab Emtansine	Breast cancer, Gastric cancer
FGFR2-Fusions	Pemigatinib, Infigratinib	Cholangiocarcinoma
FGFR3 mutation/Fusion	Erdafitinib	Urothelial Cancer
KRAS	Cetuximab, Panitumumab, Sotorasib	CRC, NSCLC
NRAS	Panitumumab	CRC
MET, exon skipping	Capmatinib	NSCLC
NTRK1, NTRK2, NTRK3-fusions	Larotrectinib, Entrectinib	Pan solid tumors
PIK3CA	Olaparib, Alpelisib	Breast cancer, Pan solid tumors
RET-Fusions	Pralsetinib, Selpercatinib	NSCLC
RET-Fusions/mutation	Selpercatinib	Thyroid cancer
Ros1-fusions	Entrectinib, Crizotinib	NSCLC

ANDERSON DIAGNOSTIC & LABS

For Enquiries, Contact:

CENTRAL PROCESSING LAB

No:38 (159), Greams Road,
Thousand Lights, Chennai-600006

+91 78248 55255

Scan here to visit
our website



UK NEQAS
Lab Reg.No.94568



info@andersondiagnostics.com